

# Did you know?

Your quest to detect hidden inflammation just got much easier

## HPT<sup>FELINE</sup>

Detect, gauge & monitor inflammatory disease

### Clinical Lab Report *example*

Feline Specific Haptoglobin (HPT)

Interpretation

HPT <sup>FELINE</sup>		REPORT DATE: 3-Sep-2014		
TEST NAME	RESULT	UNITS	FLAG	REFERENCE INTERVAL
f-HPT feline specific haptoglobin	48.9	mg/dL		Decreased (D): < 25 Optimum: 25 - 64.9 Normal: 25 - 110 Low Inflamm (L): 110.1 - 140 Mod Inflamm (M): 140.1 - 300 High Inflamm (H): ≥ 300.1

  

0	50	100	150	200	250	300	350	400	450	> 500
Dec			Low		Moderate				High	

**Flag Interpretation**  
Decreased: intravascular hemolysis, severe liver disease  
Low to Moderate: inflammatory diseases such as cancer, IBD, heart, kidney disease  
High: diseases that cause severe inflammatory response such as cancer, IBD, infection, and autoimmune conditions such as IMHA. Infection and autoimmune disorders such as IMHA often trigger HPT levels in excess of 400 mg/dL.

**General Comments**  
HPT<sup>feline</sup> measures feline-specific haptoglobin (f-HPT), an acute phase protein, produced mainly in the liver as a response to inflammation and the release of cytokines. Serum HPT has been shown to be an effective measure of general inflammation. The concentration of HPT correlates to both the severity and duration of the inflammatory stimuli.

Acute phase proteins constitute a group of proteins (e.g., CRP, haptoglobin, serum amyloid A) which are part of the innate host defense system. Their blood concentration changes rapidly in response to any tissue damaging causes such as infection, immune-mediated disorders, neoplasia, trauma, and others. HPT decreases rapidly when inflammatory stimulation is no longer present.

A low inflammatory state is healthy, however HPT has a role in removing free hemoglobin to prevent kidney damage and to return iron stores. Therefore, HPT levels can be decreased in severe liver disease and hemolytic anemia, particularly intravascular.

HPT measurements may be useful for:  
**Preventative Care Health Screen:** The high sensitivity and low specificity of the acute phase response allows for potential detection of subclinical or pre-clinical inflammatory disorders. An elevated HPT indicates the presence, but not the cause, of an inflammatory disease state.  
**Perioperative Monitoring:** HPT levels are expected to increase post-operatively, with the magnitude of increase related to the intensity of the surgical trauma. Normalization of HPT within 1-2 weeks post-operatively has been suggested as being more effective in monitoring post-surgical progress than WBC counts.  
**Monitoring Disease Progression/Response to Treatment:** HPT has been shown to be an effective marker for disease activity and response to therapy in a wide range of inflammatory conditions such as infection, inflammatory bowel disease, IMHA, heart disease, kidney disease and other systemic disorders.  
**Hemolytic Anemia:** HPT binds free hemoglobin levels will be decreased in intravascular hemolysis. In evaluating patients with suspected hemolytic anemia, HPT can be useful in differentiating intra vs extra vascular hemolysis.

Tech: RR

\*Feline haptoglobin is also available as part of the Tkfeline cancer panel. The report shown here is for feline haptoglobin as a standalone test.

 **VDI** Your Specialty Reference Lab

www.VDILab.com | (805) 577-6742

# Feline Inflammatory Disease

## Clinical Background

Chronic inflammation is both associated with, and the root-cause of disease. Left unchecked, chronic inflammation can lead to disease advancement and even death.

**HPTFELINE** is a blood test that detects systemic inflammation in response to infection, tissue damage, neoplasia, and trauma. The test is used to confirm the presence of inflammatory disease, detect hidden disorders, monitor the effectiveness of treatment, and assure complication-free recovery.

## Test Information:

Sample Type: Serum >1mL

Interferences: Gross lipemia  
Gross hemolysis

Stability: 4°C: 4 hours  
-20°C: 1 month

Reference Ranges:

**f-HPT (mg/dL):**

Decreased: <25  
Optimum: 25 – 65  
Normal: 25 – 140  
Elevated: >140

## Interpretative Information:

↑ Cancer  
Heart Disease  
Kidney Disease  
Autoimmune diseases  
Infection  
Other Inflammatory Disease

↓ Liver Disease  
Intravascular Hemolysis



## Role of inflammation

### Acute defends, chronic kills



#### Acute Inflammation

Acute (short-term) inflammation is a vital life-sustaining function. The cascade of events that occurs is needed to initiate a defense against invading organisms and to repair tissue damage that occurs from trauma, infection, and disease.

#### When inflammation becomes chronic

While acute inflammation is normally tightly controlled and part of the healing process, chronic (long-term) inflammation leads to the production of free radicals and other destructive agents, creating an environment that propagates disease.



#### Effective Wellness Screen

Cats are notorious for hiding any signs of illness. Since HPT measures inflammation, regardless of the source, the cat cannot be considered “well” if inflammation is present. Further if the inflammation becomes chronic it can lead to more serious consequences.

#### Pre-surgical test

Surgery causes tissue damage and will induce inflammation. If the cat is already in an inflammatory state further inflammation stimuli caused by surgery could induce systemic shock (SIRS). If surgery is elective, the surgery should be postponed and the cause of the inflammation investigated and resolved.

## Therapeutic Monitoring

### Track patients' response to treatment or disease progression

#### Response to Therapy

There is a direct correlation between HPT concentration and the severity and duration of the inflammation. Hence if therapy is successful, HPT will drop quickly. This is useful to monitor the treatment of infections, IBD, LSA, pancreatitis, and other diseases that induce an inflammatory response. HPT concentrations will react to changing states within the cat within hours making it a very useful measurement for “real-time” monitoring.

